


## Product datasheet for **GA207686**

### Mouse **Cmtm3** activation kit by **CRISPRa**

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | Cmtm3 CRISPRa kit - CRISPR gene activation of mouse CKLF-like MARVEL transmembrane domain containing 3   |
| Symbol:       | Cmtm3  |
| Components:   | GA207686G1, Cmtm3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA207686G2, Cmtm3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA207686G3, Cmtm3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q99LJ5</a>   |
| Synonyms:     | 9430096L06Rik; AI413895; BNAS2; Cklfsf3  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)  |

#### Product images:



[View online >](#)