


Product datasheet for **GA207691**

Mouse Upf3b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Upf3b CRISPRa kit - CRISPR gene activation of mouse UPF3 regulator of nonsense transcripts homolog B (yeast)
Symbol:	Upf3b
Components:	GA207691G1, Upf3b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207691G2, Upf3b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207691G3, Upf3b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	5730594O13Rik; AI317193; AW541158; RENT3B; UPF3X
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)