


Product datasheet for **GA207735**

Mouse Spindoc activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Spindoc CRISPRa kit - CRISPR gene activation of mouse spindlin interactor and repressor of chromatin binding
Symbol:	Spindoc
Components:	GA207735G1, Spindoc gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207735G2, Spindoc gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207735G3, Spindoc gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q05AH6</u>
Synonyms:	1700105P06Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)