


Product datasheet for **GA207740**

Mouse Mturn activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mturn CRISPRa kit - CRISPR gene activation of mouse maturin, neural progenitor differentiation regulator homolog (Xenopus)
Symbol:	Mturn
Components:	GA207740G1, Mturn gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207740G2, Mturn gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207740G3, Mturn gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8CGA4
Synonyms:	2410066E13Rik; AI448550; B230212L03Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)