

Product datasheet for **GA207764**

Mouse 9530077C05Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	9530077C05Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 9530077C05 gene
Symbol:	9530077C05Rik
Components:	GA207764G1, 9530077C05Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207764G2, 9530077C05Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207764G3, 9530077C05Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7TQE7</u>
Synonyms:	1110003N12Rik; Kiaa0895; mKIAA0895
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)