

Product datasheet for GA207772

Mouse Poglut3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Poglut3 CRISPRa kit - CRISPR gene activation of mouse protein O-glucosyltransferase 3
Symbol:	Poglut3
Components:	GA207772G1, Poglut3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207772G2, Poglut3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207772G3, Poglut3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>G5E897</u>
Synonyms:	2010004J24Rik; 4833410J10Rik; AW549401
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

