


Product datasheet for **GA207779**

Mouse Tsr3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tsr3 CRISPRa kit - CRISPR gene activation of mouse TSR3 20S rRNA accumulation
Symbol:	Tsr3
Components:	GA207779G1, Tsr3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207779G2, Tsr3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207779G3, Tsr3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5HZH2</u>
Synonyms:	0610007P22Rik; 1110014B11Rik; AL022759
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)