

Product datasheet for **GA207805**

Mouse Syncn activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Syncn CRISPRa kit - CRISPR gene activation of mouse syncollin
Symbol:	Syncn
Components:	GA207805G1, Syncn gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207805G2, Syncn gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207805G3, Syncn gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8VCK7</u>
Synonyms:	0910001K16Rik; 1810038B08Rik; SYL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)