

Product datasheet for **GA207813**

Mouse Dusp23 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dusp23 CRISPRa kit - CRISPR gene activation of mouse dual specificity phosphatase 23
Symbol:	Dusp23
UniProt ID:	<u>Q6NT99</u>
Synonyms:	1300005N15Rik; LDP-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)