

Product datasheet for **GA207843**

Mouse Wipf2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wipf2 CRISPRa kit - CRISPR gene activation of mouse WAS/WASL interacting protein family, member 2
Symbol:	Wipf2
Components:	GA207843G1, Wipf2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207843G2, Wipf2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207843G3, Wipf2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6PEV3</u>
Synonyms:	1110014J05Rik; 5730509C05Rik; AA407487; Gm1176; Wich; Wire
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)