

Product datasheet for **GA207863**

Mouse Tmed10 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmed10 CRISPRa kit - CRISPR gene activation of mouse transmembrane p24 trafficking protein 10
Symbol:	Tmed10
UniProt ID:	<u>Q9D1D4</u>
Synonyms:	1110014C03Rik; p23; p24delta1; Tmp21
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)