

Product datasheet for **GA207890**

Mouse Trpm4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Trpm4 CRISPRa kit - CRISPR gene activation of mouse transient receptor potential cation channel, subfamily M, member 4
Symbol:	Trpm4
UniProt ID:	Q7TN37
Synonyms:	1110030C19Rik; AW047689; LTrpC-4; LTRPC4; TRPM4B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)