

## Product datasheet for **GA207896**

### Mouse KansIII activation kit by CRISPRa

#### Product data:

<b>Product Type:</b>	CRISPRa Kits
<b>Description:</b>	KansIII CRISPRa kit - CRISPR gene activation of mouse KAT8 regulatory NSL complex subunit 1-like
<b>Symbol:</b>	KansIII
<b>Components:</b>	GA207896G1, KansIII gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207896G2, KansIII gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207896G3, KansIII gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
<b>UniProt ID:</b>	<u><a href="#">Q5DTI6</a></u>
<b>Synonyms:</b>	1110028C15Rik; 1190030G24; A1427598; C430010P07Rik; mKIAA4189
<b>Format:</b>	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:

