

## Product datasheet for **GA207899**

### Mouse Lce1e activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | Lce1e CRISPRa kit - CRISPR gene activation of mouse late cornified envelope 1E                                                                                                                                                                                                                             |
| Symbol:       | Lce1e                                                                                                                                                                                                                                                                                                      |
| Components:   | GA207899G1, Lce1e gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA207899G2, Lce1e gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA207899G3, Lce1e gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| Synonyms:     | 1110031B11Rik; AI507275                                                                                                                                                                                                                                                                                    |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)