

Product datasheet for **GA207911**

Mouse Ppp4r3a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ppp4r3a CRISPRa kit - CRISPR gene activation of mouse protein phosphatase 4 regulatory subunit 3A
Symbol:	Ppp4r3a
UniProt ID:	<u>Q6P2K6</u>
Synonyms:	1110034C04Rik; BC064465; mKIAA2010; Smek1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)