

Product datasheet for **GA207929**

Mouse Ms4a6d activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ms4a6d CRISPRa kit - CRISPR gene activation of mouse membrane-spanning 4-domains, subfamily A, member 6D
Symbol:	Ms4a6d
Components:	GA207929G1, Ms4a6d gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207929G2, Ms4a6d gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207929G3, Ms4a6d gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q99N07</u>
Synonyms:	1110058E16Rik; AI447446; Ms4a11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)