

Product datasheet for **GA207939**

Mouse Tmem214 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem214 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 214
Symbol:	Tmem214
Components:	GA207939G1, Tmem214 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207939G2, Tmem214 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207939G3, Tmem214 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BM55</u>
Synonyms:	1110039B18Rik; 4921530J21Rik; AI841294; AU018742; C76280
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)