


## Product datasheet for **GA207943**

### Mouse Mypn activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| Description:  | Mypn CRISPRa kit - CRISPR gene activation of mouse myopalladin                                                                                                                                                                                                                                          |
| Symbol:       | Mypn                                                                                                                                                                                                                                                                                                    |
| Components:   | GA207943G1, Mypn gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA207943G2, Mypn gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA207943G3, Mypn gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q5DTJ9</a>                                                                                                                                                                                                                                                                                  |
| Synonyms:     | 1110056A04Rik; 9830112J03; AI853556; mKIAA4170                                                                                                                                                                                                                                                          |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:



[View online »](#)