

Product datasheet for **GA207949**

Mouse Sync activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sync CRISPRa kit - CRISPR gene activation of mouse syncollin
Symbol:	Sync
UniProt ID:	<u>Q9EPM5</u>
Synonyms:	1110057H03Rik; SNIP4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)