

Product datasheet for **GA207954**

Mouse Tulp4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tulp4 CRISPRa kit - CRISPR gene activation of mouse tubby like protein 4
Symbol:	Tulp4
Components:	GA207954G1, Tulp4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207954G2, Tulp4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207954G3, Tulp4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9JIL5
Synonyms:	1110057P05Rik; 2210038L17Rik; Tusp
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)