

## Product datasheet for **GA207957**

### Mouse *Lrrn4cl* activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Lrrn4cl CRISPRa kit - CRISPR gene activation of mouse LRRN4 C-terminal like
Symbol:	Lrrn4cl
Components:	GA207957G1, Lrrn4cl gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207957G2, Lrrn4cl gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207957G3, Lrrn4cl gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3TYX2</u>
Synonyms:	1110067112Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)