

Product datasheet for **GA207973**

Mouse Zfp467 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zfp467 CRISPRa kit - CRISPR gene activation of mouse zinc finger protein 467
Symbol:	Zfp467
Components:	GA207973G1, Zfp467 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207973G2, Zfp467 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207973G3, Zfp467 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8JZL0
Synonyms:	1190001I08Rik; AB041616; AI426952; AV017623; EZI; MNCb-3350; Znf467
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)