

Product datasheet for **GA207989**

Mouse Tlcd3b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Fam57b CRISPRa kit - CRISPR gene activation of mouse family with sequence similarity 57, member B
Symbol:	Tlcd3b
Components:	GA207989G1, Tlcd3b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA207989G2, Tlcd3b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA207989G3, Tlcd3b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7TNV1</u>
Synonyms:	1500016O10Rik; A330104J06Rik; AI413816; AW060769
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)