

Product datasheet for **GA208001**

Mouse Tatdn3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tatdn3 CRISPRa kit - CRISPR gene activation of mouse TatD DNase domain containing 3
Symbol:	Tatdn3
Components:	GA208001G1, Tatdn3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208001G2, Tatdn3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208001G3, Tatdn3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3U1C6</u>
Synonyms:	1500010M24Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)