

## Product datasheet for **GA208029**

### Mouse Atp2c2 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | Atp2c2 CRISPRa kit - CRISPR gene activation of mouse ATPase, Ca++ transporting, type 2C, member 2                                                                                                                                                                                                             |
| Symbol:       | Atp2c2                                                                                                                                                                                                                                                                                                        |
| Components:   | GA208029G1, Atp2c2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA208029G2, Atp2c2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA208029G3, Atp2c2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">A7L9Z8</a>                                                                                                                                                                                                                                                                                        |
| Synonyms:     | 1810010G06Rik; mKIAA0703; Spca2                                                                                                                                                                                                                                                                               |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)