

Product datasheet for **GA208030**

Mouse Slc30a5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc30a5 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 30 (zinc transporter), member 5
Symbol:	Slc30a5
Components:	GA208030G1, Slc30a5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208030G2, Slc30a5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208030G3, Slc30a5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8R4H9
Synonyms:	1810010K08Rik; AF233321; ZnT-5; ZNT5; Zntl1; ZTL1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)