

## Product datasheet for **GA208037**

### Mouse Pnlip activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Pnlip CRISPRa kit - CRISPR gene activation of mouse pancreatic lipase
Symbol:	Pnlip
UniProt ID:	<u>Q6P8U6</u>
Synonyms:	1810007A24Rik; PL; PTL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)