

Product datasheet for **GA208044**

Mouse Kdf1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Kdf1 CRISPRa kit - CRISPR gene activation of mouse keratinocyte differentiation factor 1
Symbol:	Kdf1
Components:	GA208044G1, Kdf1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208044G2, Kdf1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208044G3, Kdf1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A2A9F4
Synonyms:	1810019J16Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)