

## Product datasheet for **GA208052**

### Mouse YdjC activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| Description:  | YdjC CRISPRa kit - CRISPR gene activation of mouse YdjC homolog (bacterial)                                                                                                                                                                                                                             |
| Symbol:       | YdjC                                                                                                                                                                                                                                                                                                    |
| Components:   | GA208052G1, YdjC gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA208052G2, YdjC gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA208052G3, YdjC gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q14BV6</a></u>                                                                                                                                                                                                                                                                           |
| Synonyms:     | 1810015A11Rik; 4930521M19Rik                                                                                                                                                                                                                                                                            |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:



[View online »](#)