

Product datasheet for **GA208073**

Mouse Sec31a activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Sec31a CRISPRa kit - CRISPR gene activation of mouse Sec31 homolog A (<i>S. cerevisiae</i>) |
| Symbol: | Sec31a |
| UniProt ID: | <u>Q3UPL0</u> |
| Synonyms: | 1810024J13Rik; ABP125; ABP130; AU042160; HSPC275; Sec3111 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)