

Product datasheet for **GA208088**

Mouse Tmem121 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem121 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 121
Symbol:	Tmem121
Components:	GA208088G1, Tmem121 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208088G2, Tmem121 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208088G3, Tmem121 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80XA0</u>
Synonyms:	2410008J05Rik; Hole
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)