

Product datasheet for **GA208094**

Mouse Plekha4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Plekha4 CRISPRa kit - CRISPR gene activation of mouse pleckstrin homology domain containing, family A (phosphoinositide binding specific) member 4
Symbol:	Plekha4
Components:	GA208094G1, Plekha4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208094G2, Plekha4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208094G3, Plekha4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	2410005C22Rik; C88169; Pepp1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)