

Product datasheet for **GA208106**

Mouse Zfp397 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zfp397 CRISPRa kit – CRISPR gene activation of mouse zinc finger protein 397
Symbol:	Zfp397
Components:	GA208106G1, Zfp397 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208106G2, Zfp397 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208106G3, Zfp397 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	2810411K16Rik; 6720480F11Rik; AI461664
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

