


Product datasheet for **GA208114**

Mouse 3300002I08Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	3300002I08Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 3300002I08 gene
Symbol:	3300002I08Rik
Components:	GA208114G1, 3300002I08Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208114G2, 3300002I08Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208114G3, 3300002I08Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	2900016H11Rik; RP23-91K11.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)