

Product datasheet for **GA208115**

Mouse Spata4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Spata4 CRISPRa kit - CRISPR gene activation of mouse spermatogenesis associated 4
Symbol:	Spata4
Components:	GA208115G1, Spata4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208115G2, Spata4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208115G3, Spata4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8K3V1</u>
Synonyms:	1700001N01Rik; Srg2; TSARG2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)