

Product datasheet for **GA208129**

Mouse Hmgb4 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Hmgb4 CRISPRa kit - CRISPR gene activation of mouse high-mobility group box 4 |
| Symbol: | Hmgb4 |
| UniProt ID: | <u>Q6P8W9</u> |
| Synonyms: | 1700001F22Rik; AI326135 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)