

Product datasheet for **GA208135**

Mouse Lelp1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lelp1 CRISPRa kit - CRISPR gene activation of mouse late cornified envelope-like proline-rich 1
Symbol:	Lelp1
Components:	GA208135G1, Lelp1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208135G2, Lelp1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208135G3, Lelp1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9DAE3</u>
Synonyms:	1700012F11Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)