

Product datasheet for **GA208145**

Mouse Slc38a4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc38a4 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 38, member 4
Symbol:	Slc38a4
Components:	GA208145G1, Slc38a4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208145G2, Slc38a4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208145G3, Slc38a4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8R1S9</u>
Synonyms:	1110012E16Rik; 1700012A18Rik; Ata3; mATA3; mNAT3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)