

Product datasheet for **GA208150**

Mouse **Cypt3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Cypt3 CRISPRa kit - CRISPR gene activation of mouse cysteine-rich perinuclear theca 3
Symbol:	Cypt3
Components:	GA208150G1, Cypt3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208150G2, Cypt3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208150G3, Cypt3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	1700006H02Rik; Ckt1r1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

