

Product datasheet for **GA208159**

Mouse Mocs3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mocs3 CRISPRa kit - CRISPR gene activation of mouse molybdenum cofactor synthesis 3
Symbol:	Mocs3
Components:	GA208159G1, Mocs3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208159G2, Mocs3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208159G3, Mocs3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A2BDX3
Synonyms:	1700020H17Rik; Uba4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)