

Product datasheet for **GA208174**

Mouse Dnajc17 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dnajc17 CRISPRa kit - CRISPR gene activation of mouse Dnaj heat shock protein family (Hsp40) member C17
Symbol:	Dnajc17
Components:	GA208174G1, Dnajc17 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208174G2, Dnajc17 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208174G3, Dnajc17 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q91WT4
Synonyms:	1700025B16Rik; C87112; D9Bwg1371e
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)