

Product datasheet for **GA208200**

Mouse Nup35 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Nup35 CRISPRa kit - CRISPR gene activation of mouse nucleoporin 35
Symbol:	Nup35
Components:	GA208200G1, Nup35 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208200G2, Nup35 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208200G3, Nup35 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8R4R6</u>
Synonyms:	35kDa; 2310006I24Rik; 5330402E05Rik; MP44; NO44
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)