

## Product datasheet for **GA208203**

### Mouse **Tsr2** activation kit by **CRISPRa**

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                         |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| <b>Description:</b>  | Tsr2 CRISPRa kit – CRISPR gene activation of mouse TSR2 20S rRNA accumulation                                                                                                                                                                                                                           |
| <b>Symbol:</b>       | Tsr2                                                                                                                                                                                                                                                                                                    |
| <b>Components:</b>   | GA208203G1, Tsr2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA208203G2, Tsr2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA208203G3, Tsr2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>UniProt ID:</b>   | <u><a href="#">Q8C8T8</a></u>                                                                                                                                                                                                                                                                           |
| <b>Synonyms:</b>     | 2310007F12Rik                                                                                                                                                                                                                                                                                           |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:

