

Product datasheet for **GA208211**

Mouse 1700030J22Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	1700030J22Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 1700030J22 gene
Symbol:	1700030J22Rik
Components:	GA208211G1, 1700030J22Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208211G2, 1700030J22Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208211G3, 1700030J22Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BHB7</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)