

Product datasheet for **GA208225**

Mouse Ripor3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ripor3 CRISPRa kit - CRISPR gene activation of mouse RIPOR family member 3
Symbol:	Ripor3
Components:	GA208225G1, Ripor3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208225G2, Ripor3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208225G3, Ripor3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A1L3T7</u>
Synonyms:	2310033K02Rik; Gm124
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)