


Product datasheet for **GA208240**

Mouse Ap5s1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ap5s1 CRISPRa kit - CRISPR gene activation of mouse adaptor-related protein 5 complex, sigma 1 subunit
Symbol:	Ap5s1
Components:	GA208240G1, Ap5s1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208240G2, Ap5s1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208240G3, Ap5s1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9D742</u>
Synonyms:	0610038L13Rik; 2310035K24Rik; AW060727; sigma5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)