

Product datasheet for GA208265

Mouse Tmbim1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmbim1 CRISPRa kit - CRISPR gene activation of mouse transmembrane BAX inhibitor motif containing 1
Symbol:	Tmbim1
Components:	GA208265G1, Tmbim1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208265G2, Tmbim1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208265G3, Tmbim1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BJZ3</u>
Synonyms:	2310061B02Rik; AA960455; AU024746; C78899; mKIAA4161; RECSI; Tmbib1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

