

Product datasheet for **GA208282**

Mouse Camsap3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Camsap3 CRISPRa kit - CRISPR gene activation of mouse calmodulin regulated spectrin-associated protein family, member 3
Symbol:	Camsap3
UniProt ID:	<u>Q80VC9</u>
Synonyms:	2310057J16Rik; Kiaa1543; Nezha
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)