

Product datasheet for **GA208283**

Mouse Slc52a3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc52a3 CRISPRa kit - CRISPR gene activation of mouse solute carrier protein family 52, member 3
Symbol:	Slc52a3
Components:	GA208283G1, Slc52a3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208283G2, Slc52a3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208283G3, Slc52a3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9D6X5
Synonyms:	2310046K01Rik; RFT2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)