

Product datasheet for **GA208288**

Mouse *Iqcg* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Iqcg</i> CRISPRa kit - CRISPR gene activation of mouse IQ motif containing G
Symbol:	<i>Iqcg</i>
Components:	GA208288G1, <i>Iqcg</i> gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA208288G2, <i>Iqcg</i> gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA208288G3, <i>Iqcg</i> gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q80W32
Synonyms:	2400003L07Rik; esgd12d; G1-374-12; repro1; stubby12d
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)